




Complete Gas Analyzer



8500-GC Gas Chromatograph

- Simultaneous analysis of He, H₂, O₂, N₂, CH₄, CO, CO₂, C₂H₄, C₂H₆, C₂H₂, H₂S from 50ppm to 100% and C₁ to C₁₂ from 1ppm to 100%
- <4 minutes up to C₆ and <10 minutes up to C₁₂
- Any sample pressure – atmospheric to 1200psi
- Moisture is no problem
- Single combined report for all analytes
- Optional software-controlled stream selection
- Optional oxygenates, alkenes, alkynes, and aromatics



Custom Analytical Solutions for Gas Analysis

Contact sales-can@scioninstruments.com for more information